



¹⁸O Labeling

Dr. Catherine Fenselau


University of Maryland, College Park, MD 20742 USA

¹⁸O₂-labeling is the Linux of isotope-labeling methods. Any laboratory can buy Water (¹⁸O, 99%) (OLM-240-99) and adapt the method to its own applications. It offers a universal strategy for uniform labeling of all peptides from any kind of protein, including modified proteins.¹ It is used to label clinical samples with unrivaled sensitivity.^{2,3} The only byproduct is water, and the immobilized catalytic enzyme can be removed mechanically. Labeling with ¹⁸O₂ is limited to binary comparisons or series thereof, and it requires a workflow with minimal manipulation of proteins since the light and heavy samples are combined at the peptide level.

Two atoms of ¹⁸O are introduced into the carboxylic acid group of every proteolytic peptide in a protein pool that has been catalyzed by members of the serine protease family, which includes trypsin, Glu-C protease, Lys-C protease and chymotrypsin. In the binding site of each protease, the residue of choice is covalently bound in

a tetrahedral intermediate, which is then disrupted by nucleophilic attack by a water molecule, cleaving the protein. The C-terminal residue in each peptide product is re-bound by the protease, e.g. Arginine and Lysine in the case of trypsin, and released by hydrolysis. If the peptide products are incubated with the catalytic enzyme in Water (¹⁸O, 99%) the level of ¹⁸O in the peptides will eventually equilibrate with the level of ¹⁸O in the solvent, preferably >95%. Peptide binding by the protease offers the advantage that cleavage of the protein can be optimized and carried out separately from labeling the peptide.⁴

Each heavy peptide weighs 4 Da more than its ¹⁶O₂ light analog. After labeling, the mixtures of heavy and light peptides are mixed, and isotope ratios of peptide pairs are determined by LC/MS. Concurrent MS/MS measurements and appropriate computer algorithms can provide peptide identification along with quantitation.



References

1. Fenselau, C.; Yao, Z. **2009**. ¹⁸O₂-Labeling in quantitative proteomic strategies: a status report. *J Proteome Res*, *8*, 2140-2143.
2. Zang, L.; Palmer Toy, D.; Hancock, W.S.; Sgroi, D.C.; Karger, B.L. **2004**. Proteomic analysis of ductal carcinoma of the breast using laser capture microdissection, LC-MS and ¹⁶O/¹⁸O isotopic labeling. *J Proteome Res*, *3*, 604-612.
3. Bantscheff, M.; Dumpelfeld, B.; Kuster, B. **2004**. Femtomol sensitivity post-digest ¹⁸O labeling for relative quantification of differential protein complex composition. *Rapid Commun Mass Spectrom*, *18*, 869-876.
4. Yao, X.; Afonso, C.; Fenselau, C. **2003**. Dissection of proteolytic ¹⁸O labeling: endoprotease-catalyzed ¹⁶O-to-¹⁸O exchange of truncated peptide substrates. *J Proteome Res*, *2*, 147-152.